

Cleaner used for Hematology Analysis

Instruction for use

PRODUCT NAME

Cleaner used for Hematology Analysis

MODEL

Cleaner

PACKING SPECIFICATIONS

View label

APPLICABLE SCOPE

It is used for hematology analyzer to perform concentrated cleaning and to be used when the machine is turned off.

TEST PRINCIPLE

The cleaning solution containing sodium hypochlorite is inhaled into the corresponding pipeline and soaked for a certain period of time to remove pipeline stains and adhering proteins, or manually clean the detection probe of the semi-automatic instrument.

MAIN COMPONENTS

Sodium hypochlorite, buffering agent.

STORAGE CONDITIONS AND EXPIRATION DATE

This product should be stored in a well ventilated room without corrosive gases and at a temperature of 2°C-38°C and a relative humidity of not more than 90%. The validity of use after opening the lid is 60 days. Storage is valid for 2 years.

APPLICABLE INSTRUMENTS

Hematology Analyzer

INSPECTION METHOD

Connect the cleaning solution to the corresponding pipeline or suck in the instrument as required by the instrument manual.

In the instrument, let the cleaning solution fill the counting pool and the corresponding pipeline.

PRODUCT PERFORMANCE INDICATORS

Appearance is light yellow, no precipitation, particles or flocculent liquid.

CAUTIONS

1. This product is a light yellow liquid, check for contamination or color change before use, and replace once found. If freezing occurs, thaw and mix well.
2. Do not take this product internally; if taken by mistake, consult a doctor promptly; avoid contact with skin or eyes; if contact with skin, rinse with water promptly; if contact with eyes, rinse with water immediately and receive medical treatment.
3. Use within the expiration date.

PRODUCTION COMPANY

Name of the applicant / manufacturer: Dezhou Guoke Medical Technology Co., Ltd

Add.: 28-31 axis,1 and 2 span, No.2 workshop, (Zhongyuan Scienceand Technology Innovation Park) No. 6596 Dongfanghong East Road,Yuanqiao Town, Economic and Technological Development Zone,Dezhou City, Shandong Province,China

E-mail: office@guoke-medical.com

European Authorized Representative:

Riomavix S.L.

Add.: Calle de Almansa 55,1D, Madrid 28039 Spain

E-mail: riomavix@gmail.com

Date of specification approval and modification date

Date of approval: 09 December 2022

PRODUCT STANDARD NO.

Product Standard No:/



【Symbols explanation】

Label	Meaning
	Date of manufacture
	In vitro diagnostic medical device
	Manufacturer
	Biological risks
	Batch code
	Temperature limit
	Use-by date
	Authorized representative in the European Community
	Volume
	Consult instructions for use
	Do not re-use
	Keep away from sunlight
	Contains sufficient for <n>tests
	Catalogue number
	CE Marking